

Technical Data Sheet

Product

**Annexin V FITC**

Cat. Number/Size

**EXB0024**

**100 tests**

**EXB0029**

**500 tests**

**For Research Use Only.**

**Not for use in diagnostic or therapeutic procedures.**

Format

FITC

Application

FC (QC tested)

Description

Annexin V FITC is intended for detection of early apoptotic cells by flow cytometry.

Contents

Annexin V FITC - 1 vial (0.5 ml / 2.5 ml) of Ready-to-use solution. No reconstitution is necessary. Reagent is provided in stabilizing phosphate buffered saline (PBS) solution containing 15mM sodium azide. Sufficient for 100 / 500 tests.

Excitation laser

blue (488 nm)

Preparation

Purified recombinant protein is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions. Unconjugated protein and free fluorochrome are removed by size-exclusion chromatography.

Formulation

Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage and handling

Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Images and References

[www.exbio.cz](http://www.exbio.cz)

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s.

EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz).

EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.