

Technical Data Sheet

Product	Anti-Hu CD344 Purified	
Cat. Number/Size	11-834-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	CD344
Clone	CH3A4A7
Format	Purified
Reactivity	Human, Non-human primates
Application	FC (QC tested)
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml
Isotype	Mouse IgG1 kappa
Specificity	The mouse monoclonal antibody CH3A4A7 recognizes an extracellular epitope of CD344 (Frizzled 4), a 7-TM protein of G-protein-coupled receptor family, which is a marker for neuronal stem cells.
Other names	Frizzled-4, FZD4, EVR1, FEVR, GPCR
Workshop	HLDA Workshop 2006
Immunogen	WERI-RB-1 retinoblastoma cells
Entrez Gene ID	8322
Gene name	FZD4
NCBI Full Gene Name	frizzled class receptor 4
UniProt ID	Q9ULV1

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.