

Technical Data Sheet

Product

Anti-Hu CD193 FITC

Cat. Number/Size

1F-161-T100

100 tests

For Research Use Only.

Not for use in diagnostic or therapeutic procedures.

Antigen	CD193
Clone	5E8
Format	FITC
Reactivity	Human
Application	FC (QC tested)
Application details	Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.
Excitation laser	blue (488 nm)
Isotype	Mouse IgG2b kappa
Specificity	The mouse monoclonal antibody 5E8 recognizes an extracellular epitope of CD193 (chemokine receptor 3), an approximately 41 kDa protein expressed above all in eosinophils and basophils.
Other names	CCR3, CKR3, CMKBR3, CC-CKR-3
Immunogen	human CD193 transfectants
Entrez Gene ID	1232
Gene name	CCR3
NCBI Full Gene Name	C-C motif chemokine receptor 3
UniProt ID	P51677

Preparation Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.

Formulation Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage and handling Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.