

Technical Data Sheet

Product **Anti-Hu CD24 PE-DyLight® 594**  
Cat. Number/Size **T5-503-T100** **100 tests**  
**For Research Use Only.**  
**Not for use in diagnostic or therapeutic procedures.**

Antigen	CD24
Clone	SN3
Format	PE-DyLight® 594
Reactivity	Human
Application	FC (QC tested)
Application details	Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 10 <sup>6</sup> cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.
Excitation laser	blue (488 nm)
Isotype	Mouse IgG1
Specificity	The antibody SN3 reacts with CD24, a 35-45 kDa heavily glycosylated cell surface antigen. CD24 is expressed by granulocytes, B lymphocytes and by some activated T cells and T cell malignancies. It is not expressed on human thymocytes.
Other names	HSA, nectadorin, SCLC4
Workshop	HLDA IV: WS Code B 136; HLDA V: WS Code B CD24.7
Immunogen	Glycoproteins purified from human NALM-1 cell line.
Entrez Gene ID	100133941
Gene name	CD24
NCBI Full Gene Name	CD24 molecule
UniProt ID	P25063

Preparation	Purified antibody is conjugated with activated tandem dye of R-phycoerythrin-DyLight®594 (PE-DyLight®594) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Formulation	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Images and References [www.exbio.cz](http://www.exbio.cz)

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.

DyLight® is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.