

Technical Data Sheet

Product	Anti-Hu SIGLEC10 Purified
Cat. Number/Size	11-932-C100 0.1 mg
	For Research Use Only.
	Not for use in diagnostic or therapeutic procedures.

Antigen	SIGLEC10
Clone	5G6
Format	Purified
Reactivity	Human
Application	IHC, ELISA, IP
Isotype	Mouse IgG1
Specificity	The mouse monoclonal antibody 5G6 recognizes an extracellular epitope of human SIGLEC10, a sialic acid-binding lectin expressed on subsets of human leucocytes.
Other names	PRO940, SGL2
Immunogen	SIGLEC10 extracellular domain fused with human IgG1 Fc fragment
Entrez Gene ID	89790
Gene name	SIGLEC10
NCBI Full Gene Name	sialic acid binding Ig like lectin 10
UniProt ID	Q96LC7

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.